

# KERN EHA 3000-1

# KERN

The compact all-round model with robust stainless steel weighing plate for use in laboratories, industries, and for teaching



Display screen size	75×32 mm
Languages of the user interface	Symbol language

## Construction

Dimension housing (W×D×H)	225×160×50 mm
Dimensions weighing platform (W×D×H)	120×120×7 mm
Dimensions weighing surface (W×D)	120×120 mm
Material housing	ABS plastic
Material display housing	Plastics
Material weighing plate	stainless steel
Material platform	Plastics
Material load cell(s)	aluminium
Levelling feet adjustable	✓

## Functions

Number of keys for operation	5
Auto-off interval(s) in battery mode/rechargeable battery mode	5 min
Tare function	manual (multi)
Counting function	✓

## Counting

Smallest piece weight when piece counting - laboratory conditions	100 mg
Smallest piece weight when piece counting - normal conditions	1 g
Counting resolution (laboratory conditions)	30.000

## Power Supply

Input voltage power supply / power [Max]	100 - 240 V
Input voltage device / power [Max]	9 V, 600 mA
Supplied plug-in power supply	Power supply & battery function (battery not included)
Plug-in power supply / adapter for countries - included with the delivery	EURO UK CH
Plug-in power supply / adapter for countries - optional	EURO UK CH
Battery	2×1.5 V AA
Battery operating time	70 h

## Environmental conditions

Humidity of environment [Max]	80 %
-------------------------------	------

## Category

Brand	KERN
Product category	Laboratory balances
Product group	School balances

## Measuring System

Readability [d]	0,1 g
Weighing capacity [Max]	3000 g
Number of weighing ranges	1
Eccentric loading at 1/3 [Max]	0,5 g
USP Minimum weight of sample (k = 2, U = 0.1%)	600 g
Reproducibility	0,3 g
Linearity	± 0,3 g
Weighing system	Strain gauge
Resolution	30.000
Tare range	3000 g
Units	g gn dwt ct oz
Default unit	g
Stabilization time	3 s
Warm-up time	120 min
Adjustment options	Adjusting with external weight
Recommended adjusting weight	1 kg (F2); 3 kg (F2)
Possible calibration points	1 kg; 3 kg

## Approval

CE mark	✓
---------	---

## Display

Display digit height	22 mm
Display backlight	✓

# KERN EHA 3000-1



The compact all-round model with robust stainless steel weighing plate for use in laboratories, industries, and for teaching

Ambient temperature [Min]	5 °C
Ambient temperature [Max]	40 °C

## Packing & Shipping

Readability force [d] (N)	1 d
Shipping method	Parcel service
Dimensions packaging (W×D×H)	245×204×76 mm
Net weight approx.	0,50 kg
Gross weight approx.	0,75 kg

## Services (optional)

Article number for adjustment at the location of installation	961-247
Article number for DAkkS calibration	963-127
Article number for certificate of conformity (verification)	969-517

## Piktogramme

### STANDARD



### OPTION

