

KERN EMB 2000-2V



School balance with integrated density determination function



Category

Brand	KERN
Product category	Laboratory balances
Product group	School balances

Measuring System

Readability [d]	0,01 g
Readability carat [d] (ct)	0,05 ct
Weighing capacity [Max]	2000 g
Weighing capacity [Max] (ct)	10000 ct
Eccentric loading at 1/3 [Max]	0,05 g
Reproducibility	0,02 g
Linearity	± 0,05 g
Weighing system	Strain gauge
Resolution	200.000
Units	g ct
Default unit	g
Stabilization time	2,5 s
Warm-up time	120 min
Adjustment options	Adjusting with external weight
Recommended adjusting weight	2 kg (F1)

Approval

CE mark	✓
---------	---

Display

Display digit height	15 mm
Display backlight	✓

Construction

Dimension housing (W×D×H)	175×250×55 mm
---------------------------	---------------

Dimensions weighing surface (Ø)	150 mm
---------------------------------	--------

Material housing	plastic
------------------	---------

Material weighing plate	plastic
-------------------------	---------

Functions

Number of keys for operation	9
Auto-off interval(s) in battery mode/rechargeable battery mode	3 min
Tare function	manual (multi)
Density determination	✓
Underfloor weighing	Hook
Interfaces	RS-232 standard

Power Supply

Input voltage power supply / power [Max]	230 V AC, 50 Hz
Input voltage device / power [Max]	9 V, 300 mA
Plug-in power supply type	Power adapter
Supplied plug-in power supply	Power supply unit
Plug-in power supply / adapter for countries - included with the delivery	EURO UK US CH
Plug-in power supply / adapter for countries - optional	US
Rechargeable battery optional	Rchrg. battery optional
Rechargeable battery	KS-A01
Rechargeable battery charging time	10 h
Rechargeable battery operating time - backlight off	50 h
Battery	9 V block
Battery operating time	12 h

Environmental conditions

Humidity of environment [Min]	5 %
Humidity of environment [Max]	80 %
Ambient temperature [Min]	5 °C
Ambient temperature [Max]	35 °C

Packing & Shipping

Readability force [d] (N)	1 d
Shipping method	Parcel service
Dimensions packaging (W×D×H)	220×290×150 mm
Net weight approx.	0,90 kg
Gross weight approx.	1,4 kg

KERN EMB 2000-2V



School balance with integrated density determination function

Services (optional)

Article number for DAkkS calibration	963-127
Article number for certificate of conformity (verification)	969-517

Piktogramme

STANDARD



OPTION

