

# KERN OBB-A 15 15



C-mount camera adapter 0,5 x



## Category

Brand	Optics
Product category	Optical instruments
Product group	Microscope camera adapter

## Construction

Camera adapter type	C-Mount
Supported camera types	Microscope Camera
C-Mount magnification	0,5 x
C-Mount focusable	✓

## Packing & Shipping

Readability force [d] (N)	1 d
Shipping method	Parcel service
Dimensions packaging (W×D×H)	100×50×50 mm
Net weight approx.	0,30 kg
Gross weight approx.	0,50 kg
Shipping weight	0,5 kg

## Piktogramme

### STANDARD

