

# KERN ODC 241



Digital microscopy brought up-to-date - Tablet with integrated camera for the very best monitoring of the sample as well as digital documentation



Plug-in power supply / adapter for countries - optional EURO AUS UK US CH

## Packing & Shipping

Readability force [d] (N)	1 d
Shipping method	Parcel service
Dimensions packaging (W×D×H)	275×230×85 mm
Net weight approx.	0,80 kg
Gross weight approx.	1,2 kg
Shipping weight	1,1 kg

## Piktogramme

### STANDARD



## Category

Brand	Optics
Product category	Microscopes
Product group	Microscope cameras

## Approval

CE mark	✓
---------	---

## Display

Display type	LCD touch
Display screen size	9,7"
Display resolution	2048×1536

## Construction

Dimension (W×D×H)	238×51×206 mm
Camera type	Tablet camera
Camera resolution	5 MP
Camera sensor type	CMOS
Camera sensor size	1/2,5"
Camera colour depth	Colour
Camera framerate	30 fps
Camera mounting types	C-mount thread

## Functions

Memory - total size	128 GB
Interfaces	WiFi standard USB 2.0 HDMI MicroSD card slot
Supported operating systems	Android

## Power Supply

Supplied plug-in power supply	Power supply unit
-------------------------------	-------------------