

KERN ODC 825

KERN

Specialists in microscopy for measurement, counting, documentation, archiving and image processing



Net weight approx.	0,60 kg
Gross weight approx.	1,00 kg
Shipping weight	1 kg

Piktogramme

STANDARD



Category

Brand	Optics
Product category	Microscopes
Product group	Microscope cameras

Approval

CE mark	✓
---------	---

Construction

Dimension (W×D×H)	80×80×45 mm
Camera type	USB2 Mikroskopkamera - C-Mount
Camera resolution	5,1 MP
Camera sensor type	CMOS
Camera sensor size	1/2,5"
Camera colour depth	Colour
Camera framerate	6,8 - 55 fps
Camera mounting types	C-mount thread

Functions

Interfaces	USB 2.0
Supported operating systems	Win XP Win Vista Win 7 Win 8 Win 10

Power Supply

Supplied plug-in power supply	no power supply, no battery
-------------------------------	--------------------------------

Packing & Shipping

Readability force [d] (N)	1 d
Shipping method	Parcel service
Dimensions packaging (W×D×H)	164×164×75 mm