

KERN ORA 32BA

KERN

Refractive index measurement for laboratories and the industry, Area of application: Sugar



Services (optional)

Article number for factory calibration	961-290
--	---------

Piktogramme

STANDARD



Category

Brand	Optics
Product category	Refractometers
Product group	Analogue refractometer

Measuring System

Measurement temperature [Min] (°C)	10 °C
Measurement temperature [Max] (°C)	30 °C
Appropriate calibration liquid	0 % (distilled water)
Scales	Brix
Measuring ranges	0 % - 32 %
Readabilites [d]	0,2 %

Construction

Dimension (W×D×H)	178×40×40 mm
-------------------	--------------

Functions

Automatic temperature compensation	✓
------------------------------------	---

Packing & Shipping

Readability force [d] (N)	1 d
Shipping method	Parcel service
Dimensions packaging (W×D×H)	215×94×65 mm
Dimensions transport container (W×D×H)	205×75×55 mm
Net weight approx.	0,30 kg
Gross weight approx.	0,35 kg
Shipping weight	0,31 kg